

# **Data Communication Protocol for MEK-8222J/K Automated Hematology Analyzer**

If you have any comments or suggestions  
on this manual, please contact us at:  
[www.nihonkohden.com](http://www.nihonkohden.com)

0614-008739A

**Copyright Notice**

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved. No part of this document may be reproduced, stored, or transmitted in any form or by any means (electronic, mechanical, photocopied, recorded, or otherwise) without the prior written permission of Nihon Kohden.

# **LICENSE AGREEMENT (Sample)**

\_\_\_\_\_ (hereafter called the “LICENSEE”) and Nihon Kohden Corporation (hereafter called LICENSOR) shall conclude an agreement reading as follows:

## **1. Definition**

The technical information subject to this agreement is “Technical Reference Manual: All communication protocol provided for connecting LICENSOR’s Product to Third party’s Product”. This technical information is to be treated as confidential knowledge between the two parties and will remain the property of LICENSEE.

## **2. Licensing**

LICENSOR shall license use of the subject technical information to LICENSEE only for the purpose of obtaining said connecting LICENSOR’s Product to Third party’s Product. LICENSEE shall not use the subject technical information for any other purpose.

## **3. Entrustment to a third party**

Regarding application and use of the subject technical information, LICENSEE shall not entrust or subcontract the subject technical information to a third party without the express written consent of LICENSOR.

## **4. Confidentiality**

LICENSEE shall maintain confidentiality of the subject technical information disclosed by LICENSOR. Also, neither LICENSEE or LICENSOR shall disclose confidential information to a third party after conclusion of this agreement.

## **5. Improved techniques**

5-1. If LICENSOR improves the subject technical information, such as through upgrading, LICENSOR shall notify LICENSEE and said improved technical information shall, after notification, become included in the subject technical information as stipulated in Article 1 of this agreement.

5-2. Any expense which becomes necessary for application of said improved technical information created by LICENSOR shall be borne by the “LICENSEE”.

## **6. Obligations**

6-1. LICENSEE shall, when using or applying the subject technical information, sufficiently learn the appropriate application of said technical information in relation to the subject equipment and take carefully planned safety measures. It is the responsibility of LICENSEE to take all necessary safety measures in regard to use of the technical information. In addition, LICENSOR shall in no way be held responsible for any death, injury or damage resulting from use of the subject technical material.

6-2. LICENSOR shall not guarantee that the use or application of the subject technical information does not infringe upon patent rights, copyrights or any other rights being owned by a third party.

6-3 When it becomes necessary to obtain permits or to make notifications in relation with use or application of the subject technical information, on the basis of the relevant national laws and regulations, LICENSEE shall obtain said permits or make said notifications at its own risk and expense.

## **7. Transfer**

LICENSEE shall not, under any circumstances, transfer any right, obligation or portion thereof arising from this agreement to a third party.

## 8. Termination

When either LICENSEE or LICENSOR infringes upon any condition of this agreement, the other party has the right to demand fulfillment of the condition or correction of the infringing act. Should the subject party fail to perform demanded fulfillment or correction within 30 days from the date of filing of said demand, the other party reserves the right to effectively terminate this agreement by written notification to the other party.

## 9. Duration

9-1. This agreement shall be valid for one year from the date of signing and concluding the present agreement unless said agreement is terminated in accordance with Article 8. The present agreement shall be automatically renewed on a yearly basis unless either party proposes otherwise.

9-2. Even after effective validity of this agreement or after authentic termination by either of the two parties to this agreement, the stipulations in Article 4, Article 6 and Paragraph 2 of the Article 8 subject circumstances cease to exist.

## 10. Procedures to be taken after discontinuation of the efficacy of this agreement

When this agreement has been canceled according to the stipulation of the aforementioned Article 8, owing to reasons attributable to the “LICENSEE”, or when the efficacy of this agreement has been discontinued, LICENSEE shall immediately discontinue use of the subject technical information and return the subject technical information disclosed by LICENSOR.

## 11. Consultations

When any discrepancy occurs in the interpretation of the stipulations of this agreement or when any discrepancy occurs between the two parties in regards to matters not converted by the stipulations of this agreement, the two parties shall discuss the matter amicably in order to reach a solution in good faith. To prove conclusion of the present agreement, two copies of this agreement shall be made and signed by each party. Each party shall maintain one copy of the signed agreements.

LICENSEE     Date:

Hospital installed:

Product Model:

Company:

(Signature):

LICENSEE's Name:

Title:

LICENSOR     Date:

Engineering Operations

Nihon Kohden Corporation,

1-31-4, Nishiochiai, Shinjuku-ku, Tokyo

(Signature):

LICENSOR's Name:

Title:

# Contents

RS-232C Data Transfer .....	1
General .....	1
Wiring (conforms to RS-232C specifications) .....	1
Transfer Format .....	2
Data Transfer Character Code .....	2
Transfer Format for PC .....	3
Distinguishing the Hematology Analyzer .....	3
Transfer Format Example (version V03-01) .....	4
Transfer Format Example (version V02-07 and V02-03) .....	12
Socket Pin Assignment .....	16
RS-232C Socket .....	16
ZK-820V Socket .....	16
USB Socket .....	16

## RS-232C Data Transfer

### General

Sample data can be transferred to the optional printer or a personal computer via the RS-232C socket on the rear panel of the hematology analyzer. One sample data is transferred each time the counting is complete. One stored sample data is transferred when the PRINT or TRANSFER key on the screen is pressed.

#### CAUTION

- In order to avoid any safety hazard, only connect personal computers which are approved by IEC 60950.
- The hematology analyzer should only be connected to an external instrument which complies with the CISPR 11 (1997), Group 1 and Class B standard.

#### CAUTION


Connect only the specified instrument to the hematology analyzer and follow the specified procedure. Failure to follow this instruction may result in electrical shock or injury to the operator, and cause fire or instrument malfunction.

#### NOTE

- To transfer the sample data to a personal computer, change the settings on the personal computer to match the data transfer conditions with this hematology analyzer beforehand.
- To transfer data to a personal computer, select PC for OUTPUT TO on the OUTPUT FORMAT screen.

#### Wiring (conforms to RS-232C specifications)

25 pin D sub connector



## Transfer Format

There are five formats for transferring data: LQ-300+, LX-300+, TM-L90, TMU295 and PC. The default setting is LQ-300+.

The following examples show the PC format.

- Data transfer starts from 02 (STX: Start of text) and ends with 03 (ETX: End of text)
- Each data separated by 0D H (CR)
- ASCII code
- Transferring order
  1. Common data block
    - Hematology analyzer information
    - Measurement data
    - Flag data
  2. Expanded data block
    - Work list data
    - Normal range setting data

### NOTE

The transfer format is fixed. Set the receiving instrument to the same format as the hematology analyzer. Before sending the data, you must prepare the receiving instrument.

### Data Transfer Character Code

Lower Bit (Hexadecimal)	Upper Bit (Hexadecimal)															
	0	1	2	3	4	5	6	7	8	9	A	B	C	D	E	F
	0		SP	0	@	P		p								
	1		!	1	A	Q	a	q								
	2		“	2	B	R	b	r								
	3		#	3	C	S	c	s								
	4		\$	4	D	T	d	t								
	5		%	5	E	U	e	u								
	6		&	6	F	V	f	v								
	7		‘	7	G	W	g	w								
	8		(	8	H	X	h	x								
	9		)	9	I	Y	i	y								
	A		*	:	J	Z	j	z								
	B		+	;	K	[	k	{								
	C		,	<	L	\	l									
	D		–	=	M	]	m	}								
	E		.	>	N	^	n	~								
	F		/	?	O	_	o									

: no character

SP: space

### Transfer Format for PC

When you set the output format to “PC” and the software version of the automated hematology analyzer is V03-01 or later, you can choose from 3 transfer formats.

- Transfer format V03-01 has flag data of SNC and PRI, expanded data, etc.
- Transfer format V02-07 is the same as V02-03. Only the sample code of the common data block is same as the V03-01.
- When the software version of the analyzer is V02-07 or before, transfer format is fixed to V02-03.

The relation of the software version and transfer format version is as follows.

The default setting is underlined.

Software version of hematology analyzer	Analysis version	Transfer format version
V01-07	V01-01	<u>V02-03</u> (Fixed)
V02-01 to V02-06	V02-01	<u>V02-03</u> (Fixed)
V02-07 to V02-09	V02-01	<u>V02-07</u> (Fixed)
V03-01 or later	V03-01	V02-03, V02-07 or <u>V03-01</u> The transfer format can be changed in the SERVICE MAINTENANCE screen. Consult your Nihon Kohden representative.

### Distinguishing the Hematology Analyzer

When connecting more than one hematology analyzer to the receiving instruments, you can distinguish the analyzer as follows.

- Transfer format V02-03 and V02-07:  
Distinguish the analyzer by the COM port No.
- Transfer format V03-01:
  - When connecting two analyzers  
You can distinguish the analyzer by the Unit No. of the extended data block. Set UNIT 1 or UNIT 2 in the SETTINGS screen. Refer to Section 3 “Changing Settings” of the operator’s manual.
  - When connecting three or more analyzers  
Distinguish the analyzer by the COM port No.



## Transfer Format Example (version V03-01)

Transfer Items	Example	No. of Byte	Transferred Data in Characters	Transferred Data in Hexadecimal
<Start> Common data block				
Start of Text	0x02	1	STX	02
Hematology analyzer information				
Type	"MEK-8222"	11	M E K - 8 2 2 2 CR	4D 45 4B 2D 38 32 32 32 20 20 0D
Parameter no.	"22"	6	2 2 CR	20 20 20 32 32 0D
Send data bytes	"01024"	6	0 1 0 2 4 CR	30 31 30 32 34 0D
Sampling mode	"CLOSED"	13	C L O S E D CR	43 4C 4F 53 45 44 20 20 20 20 20 0D
Parameter	"CBC + Diff"	13	C B C + D i f f CR	43 42 43 20 2B 20 44 69 66 66 20 0D
Sample code	"01"	3	0 1 CR	30 31 0D
Sample label	"GROUP1"	17	G R O U P 1 CR	47 52 4F 55 50 31 20 20 20 20 20 20 0D
Rack location	"01"	5	0 1 CR	20 30 31 20 0D
Seq#	"0000001"	11	0 0 0 0 0 0 1 CR	30 30 30 30 30 31 20 20 0D
Software version	"V03-03"	9	V 0 3 - 0 3 CR	56 30 33 2D 30 33 20 20 0D
Analysis program version	"V03-02"	9	V 0 3 - 0 2 CR	56 30 33 2D 30 32 20 20 0D
Format version	"V03-01"	9	V 0 3 - 0 1 CR	56 30 33 2D 30 31 20 20 0D
Total data bytes	"01536"	6	0 1 5 3 6 CR	30 31 35 33 36 0D
Data block pattern	"1"	6	1 CR	31 20 20 20 0D
Reserve data		4	SP x 3 Byte + CR x 1 Byte	20 x 3 Byte + 0D x 1 Byte
Measurement data				
Date	"20050101"	17	2 0 0 5 CR 0 1 CR 0 1 CR	CR 32 30 30 35 0D 30 31 0D 30 31 0D 20 20 20 0D
Time	"153000"	9	1 5 CR 3 0 CR 0 0 CR	31 35 0D 33 30 0D 30 30 0D
ID	"ABCDEFGH:0001"	16	A B C D E F G H : 0 0 0 1 CR	CR 41 42 43 44 45 46 47 48 3A 30 30 30 31 20 20 0D
WBC	"6.2"	7	6 . 2 CR	20 36 2E 32 20 20 0D
NE%	"70.6"	7	7 0 . 6 CR	37 30 2E 36 20 20 0D
LY%	"21.2"	7	2 1 . 2 CR	32 31 2E 32 20 20 0D
MO%	"2.5"	7	2 . 5 CR	20 32 2E 35 20 20 0D
EO%	"5.4H"	7	5 . 4 H CR	20 35 2E 34 48 20 0D
BA%	"0.3"	7	0 . 3 CR	20 30 2E 33 20 20 0D
NE	"4.4"	7	4 . 4 CR	20 34 2E 34 20 20 0D
LY	"1.3"	7	1 . 3 CR	20 31 2E 33 20 20 0D
MO	"0.2"	7	0 . 2 CR	20 30 2E 32 20 20 0D
EO	"0.2"	7	0 . 2 CR	20 30 2E 32 20 20 0D
BA	"0.0"	7	0 . 0 CR	20 30 2E 30 20 20 0D
RBC	"5.10"	7	5 . 1 0 CR	35 2E 31 30 20 20 0D
HGB	"14.4"	7	1 4 . 4 CR	31 34 2E 34 20 20 0D
HCT	"42.3"	7	4 2 . 3 CR	34 32 2E 33 20 20 0D
MCV	"86.2"	7	8 6 . 2 CR	38 36 2E 32 20 20 0D
MCH	"28.5"	7	2 8 . 5 CR	32 38 2E 35 20 20 0D
MCHC	"33.1"	7	3 3 . 1 CR	33 33 2E 31 20 20 0D
RDW	"11.5"	7	1 1 . 5 CR	31 31 2E 35 20 20 0D
PLT	"280"	7	2 8 0 CR	20 32 38 30 20 20 0D
PCT	"0.15"	7	0 . 1 5 CR	30 2E 31 35 20 20 0D
MPV	"7.2"	7	7 . 2 CR	20 37 2E 32 20 20 0D
PDW	"18.5"	7	1 8 . 5 CR	31 38 2E 35 20 20 0D
Reserve data		210	SP x 209 Byte + CR x 1 Byte	20 x 210 Byte
Flag data (with flag +, without flag (space)) (ex 20 0D)				
Leukocytosis	"+"	2	+ CR	2B 0D
Leukopenia	"+"	2	+ CR	2B 0D
Neutrophilia	"+"	2	+ CR	2B 0D
Neutropenia	"+"	2	+ CR	2B 0D
Lymphocytosis	"+"	2	+ CR	2B 0D
Lymphopenia	"+"	2	+ CR	2B 0D
Monocytosis	"+"	2	+ CR	2B 0D
Eosinophilia	"+"	2	+ CR	2B 0D
Basophilia	"+"	2	+ CR	2B 0D
Blasts	"+"	2	+ CR	2B 0D
Immature granulocyte	"+"	2	+ CR	2B 0D
Left Shift	"+"	2	+ CR	2B 0D
Atypical lymphocytes	"+"	2	+ CR	2B 0D
Poor hemolization	"+"	2	+ CR	2B 0D
Small nucleated cell	"+"	2	+ CR	2B 0D
Ly-Mo interference	"+"	2	+ CR	2B 0D
Ne-Eo interference	"+"	2	+ CR	2B 0D
Reserve data		14	SP x 13 Byte + CR x 1 Byte	20 x 13 Byte + 0D x 1 Byte
Erythrocytosis	"+"	2	+ CR	2B 0D
Anemia	"+"	2	+ CR	2B 0D
Anisocytosis	"+"	2	+ CR	2B 0D
Microcytosis	"+"	2	+ CR	2B 0D
Macrocytosis	"+"	2	+ CR	2B 0D
Hypochromia	"+"	2	+ CR	2B 0D
Abnormal MCHC	"+"	2	+ CR	2B 0D
Reserve data		10	SP x 9 Byte + CR x 1 Byte	20 x 9 Byte + 0D x 1 Byte
Thrombocytosis	"+"	2	+ CR	2B 0D
Thrombocytopenia	"+"	2	+ CR	2B 0D
PLT Clumps	"+"	2	+ CR	2B 0D
PLT-RBC interference	"+"	2	+ CR	2B 0D
Reserve data		8	SP x 7 Byte + CR x 1 Byte	20 x 7 Byte + 0D x 1 Byte
Reserve data				
Reserve data		400	SP x 399 Byte + CR x 1 Byte	20 x 399 Byte + 0D x 1 Byte
<End> Common data block				
End of Text	0x03	1	ETX	03

Transfer Items	Example	No. of Byte	Transferred Data in Characters	Transferred Data in Hexadecimal
<Start> Extended data block				
Start of Text	0x02	1	STX	02
Hematology analyzer information				
Identifier	"EXP"	4	E X P CR	45 58 50 0D
Send data bytes	"00512"	6	0 0 5 1 2 CR	30 30 35 31 32 0D
Type	"MEK-8222"	11	M E K - 8 2 2 2 CR	4D 45 4B 2D 38 32 32 32 20 20 0D
Unit no	"1"	3	1 CR	20 31 0D
Work list data				
Name	"DAVID"	27	D A V I D ... CR	44 41 56 49 44 ... 20 20 20 20 ... 0D
Sex	"MALE"	7	M A L E CR	4D 41 4C 45 20 20 0D
Date of birth	"19800219"	11	1 9 8 0 CR 0 2 CR 1 9 CR	31 39 38 30 0D 30 32 0D 31 39 0D
Age	"22"	4	2 2 CR	20 32 32 0D
Department	"INTERNAL"	14	I N T E R N A L CR	49 4E 54 45 52 4E 41 4C 20 20 20 20 0D
Physician	"WATSON"	27	W A T S O N ... CR	57 57 41 54 53 4F ... 20 20 20 ... 0D
Operator name	"STEVE"	9	S T E V E CR	53 54 45 56 45 20 20 20 0D
Comments	"No problem."	129	N o p r o b l e m . ... CR	4E 4E 6F 20 70 72 6F 62 6C 65 6D ... 20 20 20 20 20 ... 0D
Normal range table no	"0"	2	0 CR	30 0D
Work list flag	"1"	2	1 CR	31 0D
Control mode flag	"0"	2	0 CR	20 0D
Reserve data for work list data		32	SP x 31 Byte + CR x 1 Byte	20 x 31 Byte + 0D x 1 Byte
Normal range setting data				
WBC-LOW	"4.0"	5	4 . 0 CR	20 34 2E 30 0D
WBC-HIGH	"9.0"	5	9 . 0 CR	20 39 2E 30 0D
NE%-LOW	"42.0"	5	4 2 . 0 CR	34 32 2E 30 0D
NE%-HIGH	"85.0"	5	8 5 . 0 CR	38 35 2E 30 0D
LY%-LOW	"11.0"	5	1 1 . 0 CR	31 31 2E 30 0D
LY%-HIGH	"49.0"	5	4 9 . 0 CR	34 39 2E 30 0D
MO%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
MO%-HIGH	"9.0"	5	9 . 0 CR	20 39 2E 30 0D
EO%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
EO%-HIGH	"3.0"	5	3 . 0 CR	20 33 2E 30 0D
BA%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
BA%-HIGH	"2.0"	5	2 . 0 CR	20 32 2E 30 0D
NE-LOW	"1.7"	5	1 . 7 CR	20 31 2E 37 0D
NE-HIGH	"7.7"	5	7 . 7 CR	20 37 2E 37 0D
LY-LOW	"0.4"	5	0 . 4 CR	20 30 2E 34 0D
LY-HIGH	"4.4"	5	4 . 4 CR	20 34 2E 34 0D
MO-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
MO-HIGH	"0.8"	5	0 . 8 CR	20 30 2E 38 0D
EO-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
EO-HIGH	"0.3"	5	0 . 3 CR	20 30 2E 33 0D
BA-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
BA-HIGH	"0.2"	5	0 . 2 CR	20 30 2E 32 0D
RBC-LOW	"3.80"	5	3 . 8 0 CR	33 2E 38 30 0D
RBC-HIGH	"5.30"	5	5 . 3 0 CR	35 2E 33 30 0D
HGB-LOW	"11.0"	5	1 1 . 0 CR	31 31 2E 30 0D
HGB-HIGH	"17.0"	5	1 7 . 0 CR	31 37 2E 30 0D
HCT-LOW	"36.0"	5	3 6 . 0 CR	33 36 2E 30 0D
HCT-HIGH	"56.0"	5	5 6 . 0 CR	35 36 2E 30 0D
MCV-LOW	"80.0"	5	8 0 . 0 CR	38 30 2E 30 0D
MCV-HIGH	"100"	5	1 0 0 CR	20 31 30 30 0D
MCH-LOW	"28.0"	5	2 8 . 0 CR	32 38 2E 30 0D
MCH-HIGH	"36.0"	5	3 6 . 0 CR	33 36 2E 30 0D
MCHC-LOW	"31.0"	5	3 1 . 0 CR	33 31 2E 30 0D
MCHC-HIGH	"37.0"	5	3 7 . 0 CR	33 37 2E 30 0D
RDW-LOW	"11.5"	5	1 1 . 5 CR	31 31 2E 35 0D
RDW-HIGH	"16.5"	5	1 6 . 5 CR	31 36 2E 35 0D
PLT-LOW	"120"	5	1 2 0 CR	20 31 32 30 0D
PLT-HIGH	"380"	5	3 8 0 CR	20 33 38 30 0D
PCT-LOW	"0.10"	5	0 . 1 0 CR	30 2E 31 30 0D
PCT-HIGH	"1.00"	5	1 . 0 0 CR	31 2E 30 30 0D
MPV-LOW	"5.0"	5	5 . 0 CR	20 35 2E 30 0D
MPV-HIGH	"10.0"	5	1 0 . 0 CR	31 30 2E 30 0D
PDW-LOW	"12.0"	5	1 2 . 0 CR	31 32 2E 30 0D
PDW-HIGH	"18.0"	5	1 8 . 0 CR	31 38 2E 30 0D
<End> Extended data block				
End of Text	0x03	1	ETX	03

- "0x0D" is at the end of each item as separating character.
- Each item is entered by ASCII code. Only "Start of Text" and "End of Text" are entered by control code.
- Each item is initialized at space (0x20). Unused items are also initialized at space and have "0x0D" in the end.

Transferring Item Description  
Common data block

Item No.	Item	No. of Bytes	Description																
Hematology analyzer information																			
1	Type	11	Sets the model of the hematology analyzer "MEK-8222".																
2	Parameter No.	6	Sets the number of measurable parameters "22".																
3	Send data bytes	6	Set the total number of transferring data "01024".																
4	Sampling mode	13	<p>Sets the sampling method.</p> <p>Closed mode: "CLOSED"</p> <p>Closed high dilution mode 1: "HIGH (C)"</p> <p>Closed high dilution mode 2: "HIGHER (C)"</p> <p>Open mode: "MANUAL"</p> <p>Open high dilution mode 1: "HIGH (M)"</p> <p>Open high dilution mode 2: "HIGHER (M)"</p> <p>Open low dilution mode: "LOW (M)"</p> <p>Pre-dilution mode (20 µL): "PreDIL 20"</p> <p>Pre-dilution mode (10 µL): "PreDIL 10"</p>																
5	Parameter	13	<p>Sets the measuring parameters.</p> <p>When measuring 22: CBC + Diff</p> <p>When measuring 8: CBC</p>																
6	Sample code	3	<p>Sets the code for sample type.</p> <p>When measuring the normal sample, the No. of the normal range is set.</p> <table border="1"> <thead> <tr> <th>Normal range name</th><th>No.</th><th>Normal range name</th><th>No.</th></tr> </thead> <tbody> <tr> <td>GROUP 1</td><td>01</td><td>GROUP 4</td><td>04</td></tr> <tr> <td>GROUP 2</td><td>02</td><td>GROUP 5</td><td>05</td></tr> <tr> <td>GROUP 3</td><td>03</td><td></td><td></td></tr> </tbody> </table> <ul style="list-style-type: none"> <li>X-R (NORMAL) or L &amp; J (NORMAL): "21" (first time) or "22" (second time)</li> <li>X-R (LOW) or L &amp; J (LOW): "23" (first time) or "24" (second time)</li> <li>X-R (HIGH) or L &amp; J (HIGH): "25" (first time) or "26" (second time)</li> </ul>	Normal range name	No.	Normal range name	No.	GROUP 1	01	GROUP 4	04	GROUP 2	02	GROUP 5	05	GROUP 3	03		
Normal range name	No.	Normal range name	No.																
GROUP 1	01	GROUP 4	04																
GROUP 2	02	GROUP 5	05																
GROUP 3	03																		
7	Sample label	17	<p>Sets the sample type.</p> <p>When measuring the normal sample, the name of the normal range (GROUP 1 to GROUP 5) is set.</p> <ul style="list-style-type: none"> <li>X-R (NORMAL) or L &amp; J (NORMAL): "X-R NORMAL"</li> <li>X-R (LOW) or L &amp; J (LOW): "X-R LOW"</li> <li>X-R (HIGH) or L &amp; J (HIGH): "X-R (HIGH)"</li> </ul>																

8	Rack location	5	<p>Sets the rack ID and rack position.</p> <p>1 Bytes: rack ID: “1” to “9” Emergency sample: “E” Open mode: “M” Not use rack ID (other than emergency sample and open mode): space (0x20)</p> <p>2 to 3 bytes: rack position: “01” to “50” Open mode: “MM”</p> <p>4 bytes: space (0x20)</p> <p>5 bytes: CR (0x0D)</p> <p>&lt;Setting example&gt;</p> <ul style="list-style-type: none"> <li>• Rack position 10: 10</li> <li>• Emergency sample, rack position 1: “E01”</li> <li>• Open mode: “MMM”</li> </ul>
9	Seq#	11	Sets the sequence number.
10	Software version	9	Sets the software version.
11	Analysis program version	9	Sets the analysis program version.
12	Format version	9	Sets the RS-232C output format version.
13	Total data bytes	6	Sets the all sent data size “01536” including common data block and extended data block.
14	Data clock pattern	6	Sets the data block that is send after common data block. When the extended data block is send, “1” is set.
15	Reserve data	4	Sets space (0x20) × 3 + CR (0x0D)..
Measurement data			
1	Date	17	Sets date. Year, month and day are separated by CR (0x0D). A day of the week is not set. Space (0x20) is set.
2	Time	9	Sets time. Hour, minute and second are separated by CR (0x0D).
3	ID	16	Sets the sample ID number.
4	WBC to PDW	7 each	<p>Sets as follows.</p> <p>1 to 4 bytes: measurement value when over the limits: “OVER” when measurement is impossible: space (0x20)</p> <p>5 to 6 bytes: abnormal value mark</p> <p>7 bytes: CR (0x0D)</p>
5	Reserve data	210	Sets space (0x20) × 209 + CR (0x0D).
Flag data			
1	Leukocytosis to Ne-Eo interference	2 each	Sets “+” when there is a flag, space (0x20) when there is no flag.
2	Reserve data	14	Sets space (0x20) × 13 + CR (0x0D).
3	Erythrocytosis to abnormal MCHC	2 each	Sets “+” when there is a flag, space (0x20) when there is no flag.
4	Reserve data	10	Sets space (0x20) × 9 + CR (0x0D).
5	Thrombocytosis to PLT-RBC interference	2 each	Sets “+” when there is a flag, space (0x20) when there is no flag.
6	Reserve data	8	Sets space (0x20) × 7 + CR (0x0D).
Reserved data			
1	Reserve data	400	Sets space (0x20) × 399 + CR (0x0D).

#### NOTE

In the reserved data, new items by character codes may be set when the software is upgraded in the future.

# Extended data block

Item No.	Item	No. of Bytes	Description
Hematology analyzer information			
1	Identifier	4	Sets the data block identifier "EXP".
2	Send data bytes	6	Sets the total number of transferred extended data block "00512".
3	Type	11	Sets the model of the hematology analyzer "MEK-8222".
4	Unit No.	3	Sets the analyzer name. When set to "UNIT 1": 1 When set to "UNIT 2": 2
Work list data			
1	Name	27	Sets the patient name. When not using the work list, a space (0x20) × 26 + CR (0x0D) is set.
2	Sex	7	Sets the patient sex, MALE or FEMALE. When not using the work list, a space (0x20) × 6 + CR (0x0D) is set.
3	Date of birth	11	Sets the patient date of birth. When not using the work list, a space (0x20) × 10 + CR (0x0D) is set.
4	Age	4	Sets the patient age. When not using the work list, a space (0x20) × 3 + CR (0x0D) is set.
5	Department	14	Sets the department of the patient. When not using the work list, a space (0x20) × 13 + CR (0x0D) is set.
6	Physician	27	Sets the name of the physician. When not using the work list, a space (0x20) × 26 + CR (0x0D) is set.
7	Operator name	9	Sets the name of the operator.
8	Comments	129	Sets the comment. When not using the work list, a space (0x20) × 128 + CR (0x0D) is set.
9	Normal range table no.	2	Sets the group number of the normal range used for measurement. GROUP 1: 0 GROUP 2: 1 GROUP 3: 2 GROUP 4: 3 GROUP 5: 4
10	Work list flag	2	Sets the work list flag. When measured using the work list: 1 When measured without work list: 0
11	Control mode flag	2	Sets space (0x20)+CR (0x0D).
12	Reserve data	32	Sets space (0x20) × 31 + CR (0x0D).
Normal range data			
1	WBC-HIGH to PDW-LOW	5 each	Sets 4 bytes + CR (0x0D) for normal range data used for measurement.

## Abnormal mark

Abnormal mark shows following error.

Item	Mark	Description	Item	Mark	Description
All	H	Over normal range	LY	?	Room temp high, room temp low or optical count error
	L	Under normal range		*	Blasts, atypical Ly or Ly-Mo interference
WBC	?	WBC sample error	LY%	?	Room temp high, room temp low or optical count error
	!	Poor hemolyzation		*	Blasts, atypical Ly or Ly-Mo interference
	C	PLT clumps	MO	?	Optical count error
	*	Small nucleated cell or WBC previous data OVER		*	Blasts, immature Gr, left shift, atypical Ly or Ly-Mo interference
RBC	?	RBC sample error	MO%	?	Optical count error
	*	RBC-PLT interference or RBC previous data OVER		*	Blasts, immature Gr, left shift, atypical Ly or Ly-Mo interference
HGB	!	HGB voltage high	EO	?	Optical count error
	?	HGB voltage low		*	Blasts or Ne-Eo interference
	*	HGB circuit error, WBC OVER or WBC/HGB previous data OVER	EO%	?	Optical count error
HCT	*	HCT previous data OVER		*	Blasts or Ne-Eo interference
MCHC	!	RBC sample error	BA	?	Optical count error
PLT	?	RBC sample error		*	Blasts, immature Gr or left shift
	C	PLT clumps	BA%	?	Optical count error
	*	RBC-PLT interference or PLT previous data OVER		*	Blasts, immature Gr or left shift
NE	?	Room temp high, room temp low or optical count error			
	*	Blasts, immature Gr, left shift or Ne-Eo interference			
NE%	?	Room temp high, room temp low or optical count error			
	*	Blasts, immature Gr, left shift or Ne-Eo interference			

## Data indication

Data indication (alarm/flag) and their description is as follows.

Classification		Description
1	Data cannot be analyzed	Data cannot be analyzed.
2	Measurement alarm (data not displayed)	Error found during measurement.
3	Measurement alarm (data displayed)	Measurement error due to surrounding temperature out of specified range. Measured data is displayed but measurement accuracy is not reliable.
4	Out of measuring range	Out of measuring range.
5	Data with low reliability	Abnormal flag detected in the sample. Measurement accuracy is not reliable due to abnormal cell. <ul style="list-style-type: none"> <li>When WBC flag appears, all WBC parameters are affected by the abnormal cell. “*” is displayed beside the parameter which is greatly affected.</li> <li>When there is possibility of PLT coagulation, “C” is displayed beside the parameter.</li> <li>When there is possibility of poor hemolyzation, “!” is displayed beside the parameter.</li> </ul>
6	Out of normal range	Out of normal range setting.

## Data, alarm, flag and data transfer

	Classification	Data display	Symbol display	Alarm message/ flag display	Data transfer
1	Data cannot be analyzed	None	None	None	None
2	Measurement alarm (data not displayed)	None	None	Alarm code no. and message	Abbreviation of the alarm
3	Measurement alarm (data displayed)	Data displayed	“?” beside numeric data	Alarm code no. and message	Data and identifier
4	Out of measuring range	“OVER” message displayed	None	None	Abbreviation of “OVER”
5	Data with low reliability	Data displayed	“*”, “!” or “C” beside numeric data	Flag name	Data and identifier
6	Out of normal range	Data displayed	“H” or “L” beside numeric data	None	Data and identifier

## Abbreviation of the alarm and example of the RS-232C data output format

Example: PLT clumps

WBC LEVEL 1

Transfer code: 0x4C, 0x45, 0x56, 0x45, 0x4C, 0x31, 0x0D

Alarm name	code	WBC	RBC
WBC fluid level 1	A021	LEVEL 1	
WBC fluid level 2	A022	LEVEL 2	
WBC fluid level 3	A023	LEVEL 3	
WBC bubble 1	A024	BBL 1	
WBC bubble 2	A025	BBL 2	
WBC bubble 3	A026	BBL 3	
WBC bubble 4	A027	BBL 4	
WBC clogged	A029	CLOG	
WBC hardware noise	A031	NOISE 2	
WBC software noise	A032	NOISE 1	
RBC fluid level 1	A041		LEVEL 1
RBC fluid level 2	A042		LEVEL 2
RBC fluid level 3	A043		LEVEL 3
RBC bubble 1	A044		BBL 1
RBC bubble 2	A045		BBL 2
RBC bubble 3	A046		BBL 3
RBC bubble 4	A047		BBL 4
RBC clogged	A049		CLOG
RBC hardware noise	A051		NOISE 2
RBC software noise	A052		NOISE 1

## Explanation of low reliable indicator and example of RS-232C data output

Example: PLT clumps

WBC 7.3C

Transfer code: 0x20, 0x37, 0x2E, 0x33, 0x43, 0x20, 0x0D

PLT 280C

Transfer code: 0x20, 0x32, 0x38, 0x30, 0x43, 0x20, 0x0D

Flag	Flag Class.	Parameters											
		NE% NE#	LY% LY#	MO% MO#	EO% EO#	BA% BA#	WBC	RBC	HGB	HCT	MCHC	PLT	PCT MPV PDW
Blasts	WBC	*	*	*	*	*							
Immature Granulocyte		*		*		*							
Left Shift		*		*		*							
Atypical Lymphocytes			*	*									
Small Nucleated Cell							*						
Ly-Mo Interference			*	*									
Ne-Eo Interference		*			*								
PLT-RBC Interference	RBC/PLT							*				*	
PLT Clumps							C					C	C
WBC previous data OVER	Other						*		*				
RBC previous data OVER								*					
HGB previous data OVER									*				
HCT previous data OVER										*			
PLT previous data OVER												*	
WBC OVER									*				
Abnormal MCHC	Specimen										!		
Poor hemolyzation							!						

Alarm message	Code No.	Parameters											
		NE% NE#	LY% LY#	MO% MO#	EO% EO#	BA% BA#	WBC	RBC	HGB	HCT	MCHC	PLT	
WBC sample error	A030						?						
RBC sample error	A050							?					?
HGB voltage low	A061								?				
HGB voltage high	A062								!				
HGB circuit error	A063								*				
Room temp high	A091	?	?										
Room temp low	A092	?	?	?									
Optical count error	A082	?	?	?	?	?							



## Transfer Format Example (version V02-07 and V02-03)

Transfer Items	Example	No. of Byte	Transferred Data in Characters	Transferred Data in Hexadecimal
<Start> Common data block				
Start of Text	0x02	1	STX	02
Hematology analyzer information				
Type	"MEK-8222"	11	M E K - 8 2 2 2 CR	4D 45 4B 2D 38 32 32 32 20 0D
Parameter no.	"22"	6	2 2 CR	20 20 20 32 32 0D
Send data bytes	"1024"	6	0 1 0 2 4 CR	30 31 30 32 34 0D
Sample type	"CLOSED"	13	C L O S E D CR	43 4C 4F 53 45 44 20 20 20 20 0D
Parameter	"CBC + Diff"	13	C B C + D i f f CR	43 42 43 20 2B 20 44 69 66 66 20 0D
Sample code	"00"	3	0 0 CR	30 30 0D
Sample type name	""	17	CR	20 20 20 20 20 20 20 20 20 20 20 20 20 20 20 20 0D
Rack location	"01"	5	0 1 CR	30 31 20 20 0D
Seq#	"0000001"	11	0 0 0 0 0 0 1 CR	30 30 30 30 30 30 31 20 20 20 0D
Reserve data		43	SP x 43 Byte	20 x 43 Byte
Measurement data				
Date	"20020725"	17	2 0 0 2 CR 0 7 CR 2 5 CR	32 30 30 32 0D 30 37 0D 32 35 0D 20 20 20 20 0D
Time	"153000"	9	1 5 CR 3 0 CR 0 0 CR	31 35 0D 33 30 0D 30 30 0D
ID	"ABCDEFGH:0001"	16	A B C D E F G H : 0 0 0 1 CR	41 42 43 44 45 46 47 48 3A 30 30 30 31 20 20 0D
WBC	"6.2"	7	6 . 2 CR	20 36 2E 32 20 20 0D
NE%	"70.6"	7	7 0 . 6 CR	37 30 2E 32 20 20 0D
LY%	"21.2"	7	2 1 . 2 CR	32 31 2E 32 20 20 0D
MO%	"2.5"	7	2 . 5 CR	20 32 2E 35 20 20 0D
EO%	"5.4H"	7	5 . 4 H CR	20 35 2E 34 48 20 0D
BA%	"0.3"	7	0 . 3 CR	20 30 2E 33 20 20 0D
NE	"4.4"	7	4 . 4 CR	20 34 2E 34 20 20 0D
LY	"1.3"	7	1 . 3 CR	20 31 2E 33 20 20 0D
MO	"0.2"	7	0 . 2 CR	20 30 2E 32 20 20 0D
EO	"0.2"	7	0 . 2 CR	20 30 2E 32 20 20 0D
BA	"0.0"	7	0 . 0 CR	20 30 2E 30 20 20 0D
RBC	"5.10"	7	5 . 1 0 CR	35 2E 31 30 20 20 0D
HGB	"14.4"	7	1 4 . 4 CR	31 34 2E 34 20 20 0D
HCT	"42.3"	7	4 2 . 3 CR	34 32 2E 33 20 20 0D
MCV	"86.2"	7	8 6 . 2 CR	38 36 2E 32 20 20 0D
MCH	"28.5"	7	2 8 . 5 CR	32 38 2E 35 20 20 0D
MCHC	"33.1"	7	3 3 . 1 CR	33 33 2E 31 20 20 0D
RDW	"11.5"	7	1 1 . 5 CR	31 31 2E 35 20 20 0D
PLT	"280"	7	2 8 0 CR	20 32 38 30 20 20 0D
PCT	"0.15"	7	0 . 1 5 CR	30 2E 31 35 20 20 0D
MPV	"7.2"	7	7 . 2 CR	20 37 2E 32 20 20 0D
PDW	"18.5"	7	1 8 . 5 CR	31 38 2E 35 20 20 0D
Reserve data		210	SP x 210 Byte	20 x 210 Byte
Flag data (with flag +, without flag (space))				
Leukocytosis	"+"	2	+ CR	2B 0D
Leukopenia	"+"	2	+ CR	2B 0D
Neutrophilia	"+"	2	+ CR	2B 0D
Neutropenia	"+"	2	+ CR	2B 0D
Lymphocytosis	"+"	2	+ CR	2B 0D
Lymphopenia	"+"	2	+ CR	2B 0D
Monocytosis	"+"	2	+ CR	2B 0D
Eosinophilia	"+"	2	+ CR	2B 0D
Basophilia	"+"	2	+ CR	2B 0D
Blasts	"+"	2	+ CR	2B 0D
Immature	"+"	2	+ CR	2B 0D
Left Shift	"+"	2	+ CR	2B 0D
Atypical	"+"	2	+ CR	2B 0D
WBC hemolysing	"+"	2	+ CR	2B 0D
Reserve data		20	SP x 20 Byte	20 x 20 Byte
Erythrocytosis	"+"	2	+ CR	2B 0D
Anemia	"+"	2	+ CR	2B 0D
Anisocytosis	"+"	2	+ CR	2B 0D
Microcytosis	"+"	2	+ CR	2B 0D
Macrocytosis	"+"	2	+ CR	2B 0D
Hypochromia	"+"	2	+ CR	2B 0D
MCHC error	"+"	2	+ CR	2B 0D
Reserve data		10	SP x 10 Byte	20 x 10 Byte
Thrombocytosis	"+"	2	+ CR	2B 0D
Thrombocytopenia	"+"	2	+ CR	2B 0D
PLT Clumps	"+"	2	+ CR	2B 0D
Reserve data		10	SP x 10 Byte	20 x 10 Byte
Reserve data				
Reserve data		400	SP x 400 Byte	20 x 400 Byte
<End> Common data block				
End of Text	0x03	1	ETX	03

Transfer Items	Example	No. of Byte	Transferred Data in Characters	Transferred Data in Hexadecimal
<Start> Extended data block				
Start of Text	0x02	1	STX	02
Hematology analyzer information				
Identifier	"EXP"	4	E X P CR	45 58 50 0D
Send data bytes	"512"	6	0 0 5 1 2 CR	30 30 35 31 32 0D
Type	"MEK-8222"	11	M E K - 8 2 2 2 CR	4D 45 4B 2D 38 32 32 32 20 20 0D
Unit no	"1"	3	1 CR	20 31 0D
Work list data				
Name	"DAVID"	27	D A V I D ... CR	44 41 56 49 44 ... 20 20 20 20 ... 0D
Sex	"MALE"	7	M A L E CR	4D 41 4C 45 20 20 0D
Date of birth	"19800219"	11	1 9 8 0 CR 0 2 CR 1 9 CR	31 39 38 30 0D 30 02 0D 31 39 0D
Age	"22"	4	2 2 CR	20 32 32 0D
Department	"INTERNAL"	14	I N T E R N A L CR	49 4E 54 45 52 4E 41 4C 20 20 20 20 0D
Physician	"WATSON"	27	W A T S O N ... CR	57 41 54 53 4F 4E ... 20 20 20 ... 0D
Operator name	"STEVE"	9	S T E V E CR	53 54 45 56 45 20 20 20 0D
Comments	"No problem."	129	N o p r o b l e m . ... CR	4E 6F 20 70 72 6F 62 6C 65 6D ... 20 20 20 20 20 ... 0D
Normal range table	"0"	2	0 CR	30 0D
Work list flag	"0"	2	0 CR	30 0D
Control mode flag	"0"	2	0 CR	20 0D
Reserve data for work list data		32	SP x 32 Byte	20 x 32 Byte
Normal range setting data				
WBC-LOW	"4.0"	5	4 . 0 CR	20 34 2E 30 0D
WBC-HIGH	"9.0"	5	9 . 0 CR	20 39 2E 30 0D
NE%-LOW	"42.0"	5	4 2 . 0 CR	34 32 2E 30 0D
NE%-HIGH	"85.0"	5	8 5 . 0 CR	38 35 2E 30 0D
LY%-LOW	"11.0"	5	1 1 . 0 CR	31 31 2E 30 0D
LY%-HIGH	"49.0"	5	4 9 . 0 CR	34 39 2E 30 0D
MO%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
MO%-HIGH	"9.0"	5	9 . 0 CR	20 39 2E 30 0D
EO%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
EO%-HIGH	"3.0"	5	3 . 0 CR	20 33 2E 30 0D
BA%-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
BA%-HIGH	"2.0"	5	2 . 0 CR	20 32 2E 30 0D
NE-LOW	"1.7"	5	1 . 7 CR	20 31 2E 37 0D
NE-HIGH	"7.7"	5	7 . 7 CR	20 37 2E 37 0D
LY-LOW	"0.4"	5	0 . 4 CR	20 30 2E 34 0D
LY-HIGH	"4.4"	5	4 . 4 CR	20 34 2E 34 0D
MO-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
MO-HIGH	"0.8"	5	0 . 8 CR	20 30 2E 38 0D
EO-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
EO-HIGH	"0.3"	5	0 . 3 CR	20 30 2E 33 0D
BA-LOW	"0.0"	5	0 . 0 CR	20 30 2E 30 0D
BA-HIGH	"0.2"	5	0 . 2 CR	20 30 2E 32 0D
RBC-LOW	"3.80"	5	3 . 8 0 CR	33 2E 38 30 0D
RBC-HIGH	"5.30"	5	5 . 3 0 CR	35 2E 33 30 0D
HGB-LOW	"11.0"	5	1 1 . 0 CR	31 31 2E 30 0D
HGB-HIGH	"17.0"	5	1 7 . 0 CR	31 37 2E 30 0D
HCT-LOW	"36.0"	5	3 6 . 0 CR	33 36 2E 30 0D
HCT-HIGH	"56.0"	5	5 6 . 0 CR	35 36 2E 30 0D
MCV-LOW	"80.0"	5	8 0 . 0 CR	38 30 2E 30 0D
MCV-HIGH	"100"	5	1 0 0 CR	20 31 30 30 0D
MCH-LOW	"28.0"	5	2 8 . 0 CR	32 38 2E 30 0D
MCH-HIGH	"36.0"	5	3 6 . 0 CR	33 36 2E 30 0D
MCHC-LOW	"31.0"	5	3 1 . 0 CR	33 31 2E 30 0D
MCHC-HIGH	"37.0"	5	3 7 . 0 CR	33 37 2E 30 0D
RDW-LOW	"11.5"	5	1 1 . 5 CR	31 31 2E 35 0D
RDW-HIGH	"16.5"	5	1 6 . 5 CR	31 36 2E 35 0D
PLT-LOW	"120"	5	1 2 0 CR	20 31 32 30 0D
PLT-HIGH	"380"	5	3 8 0 CR	20 33 38 30 0D
PCT-LOW	"0.10"	5	0 . 1 0 CR	30 2E 31 30 0D
PCT-HIGH	"1.00"	5	1 . 0 0 CR	31 2E 30 30 0D
MPV-LOW	"5.0"	5	5 . 0 CR	20 35 2E 30 0D
MPV-HIGH	"10.0"	5	1 0 . 0 CR	31 30 2E 30 0D
PDW-LOW	"12.0"	5	1 2 . 0 CR	31 32 2E 30 0D
PDW-HIGH	"18.0"	5	1 8 . 0 CR	31 38 2E 30 0D
<End> Common data block				
End of Text	0x03	1	ETX	03

## Transferring Item Description

## Common data block

Item No.	Item	No. of Bytes	Description
Hematology analyzer information			
1	Type	11	Sets the model of the hematology analyzer "MEK-8222".
2	Parameter no	6	Sets the number of measurable parameters "22".
3	Send data bytes	6	Set the total number of transferring data "1024".
4	Sample type	13	Sets the sampling method. Closed mode: "CLOSED" Manual mode: "MANUAL" Capillary mode: "CAPILLARY" High WBC mode: "HIGH WBC" Low WBC mode: "LOW WBC"
5	Parameter	13	Sets the measuring parameters. When measuring 22: CBC + Diff When measuring 8: CBC
6	Sample code	3	Sets the code for sample type. When measuring hematology control: "00"
7	Sample type name	17	Sets the sample type.
8	Rack location	5	Sets the sample rack location number.
9	Seq#	11	Sets the sequence number.
10	Reserve data	43	Sets space (0x20)×43.
Measurement data			
1	Date	17	Sets date. Year, month and day are separated by CR (0x0D).
2	Time	9	Sets time. Hour, minute and second are separated by CR (0x0D).
3	ID	16	Sets the sample ID number.
4	WBC to PDW	7 each	Sets measured data (4 bytes) and flag (2 bytes).
5	Reserve data	210	Sets space (0x20)×210.
Flag data			
1	Leukocytosis to Atypical	2 each	Sets "+" when there is a flag, space (0x20) when there is no flag.
2	Reserve data	20	Sets space (0x20)×20.
3	Erythrocytosis to MCHC error	2 each	Sets "+" when there is a flag, space (0x20) when there is no flag.
4	Reserve data	10	Sets space (0x20)×10.
5	Thrombocytosis to PLT Clumps	2 each	Sets "+" when there is a flag, space (0x20) when there is no flag.
6	Reserve data	10	Sets space (0x20)×10.
Reserved data			
1	Reserve data	400	Sets space (0x20)×400.

**NOTE**

- It is recommended to program the PC to receive data from the hematology analyzer so that the received items are differentiated by the No. of Bytes for each item and the total no. of bytes and not by CR.
- New sample codes may be added when the software is upgraded.
- In the reserved data, new items by character codes may be set when the software is upgraded in the future.

# Extended data block

Item No.	Item	No. of Bytes	Description
Hematology analyzer information			
1	Identifier	4	Sets the data block identifier "EXP".
2	Send data bytes	6	Sets the total number of transferred extended data block "512".
3	Type	11	Sets the model of the hematology analyzer "MEK-8222".
4	Unit No.	3	Sets the analyzer name. When set to "UNIT 1": 1 When set to "UNIT 2": 2
Work list data			
1	Name	27	Sets the patient name. When not using the work list, a space (0x20) is set.
2	Sex	7	Sets the patient sex, MALE or FEMALE. When not using the work list, a space (0x20) is set.
3	Date of birth	11	Sets the patient date of birth. When not using the work list, a space (0x20) is set.
4	Age	4	Sets the patient age. When not using the work list, a space (0x20) is set.
5	Department	14	Sets the department of the patient. When not using the work list, a space (0x20) is set.
6	Physician	27	Sets the name of the physician. When not using the work list, a space (0x20) is set.
7	Operator name	9	Sets the name of the operator.
8	Comments	129	Sets the comment. When not using the work list, a space (0x20) is set.
9	Normal range table no.	2	Sets the group number of the normal range used for measurement. GROUP 1: 0 GROUP 2: 1 GROUP 3: 2 GROUP 4: 3 GROUP 5: 4
10	Work list flag	2	Sets the work list flag. When measured using the work list: 1 When measured without work list: 0
11	Control mode flag	2	Sets space (0x20)+CR (0x0D).
12	Reserve data	32	Sets space (0x20)×32.
Normal range data			
1	WBC-HIGH to PDW-LOW	5 each	Sets 4 bytes for normal range data used for measurement.

## Socket Pin Assignment


### CAUTION

Connect only the specified instrument to the hematology analyzer and follow the specified procedure. Failure to follow this instruction may result in electrical shock or injury to the operator, and cause fire or instrument malfunction.

### RS-232C Socket

D sub 25 pin (female)

SDBB-25S-SL-LNK (02)



Pin No.	Signal	Pin No.	Signal	Pin No.	Signal
1	NC	10	NC	19	NC
2	TXD	11	NC	20	(To pin 6)
3	RXD	12	NC	21	NC
4	RTS	13	NC	22	NC
5	CTS	14	NC	23	NC
6	(To pin 20)	15	NC	24	NC
7	SG	16	NC	25	NC
8	NC	17	NC		
9	NC	18	NC		

### ZK-820V Socket



Pin No.	Signal	Pin No.	Signal	Pin No.	Signal
1	FG	4	(To pin 6)	7	CTS
2	TXD	5	SG	8	RTS
3	RXD	6	(To pin 4)	9	VCC

### USB Socket



Pin No.	Signal
1	NC
2	USB-D+
3	USB-D-
4	ED